

# BIOLIFT

SKIN  
BOOSTER



# Biolife-Skinbooster

## The New Skin Rejuvenation Revolution

Biolife-Skinbooster is a groundbreaking facial renewal solution based on stabilised high-molecular-weight hyaluronic acid. This innovative treatment can completely remodel the skin, restore moisture balance, improve skin structure and elasticity, and stimulate collagen regeneration. Through micro-injections of a gel-like hyaluronic acid, Biolife-Skinbooster can significantly enhance skin tightness, elasticity and radiance, not only preventing skin ageing but also repairing skin damaged by excessive sun exposure, revealing a rejuvenated youthful complexion.

# Unique Technology of Biolife-Skinbooster

Biolife-Skinbooster uses a unique and stable high-molecular-weight hyaluronic acid technology that can deeply penetrate the skin, providing long-lasting hydration, improving elasticity and radiance. The initial treatment stage adopts micro-injection, gradually restoring the skin's moisture balance and collagen production, effectively reducing wrinkles and pigmentation, and delaying skin ageing.

## Biolife-Skinbooster

Unlike traditional hyaluronic acid injection treatments, Biolife-Skinbooster uses a patented formulation of soft gel-like hyaluronic acid that can slowly release active ingredients in the skin, continuously improving elasticity and radiance. The unique micro-injection method can minimise irritation, allowing the skin to gradually adapt and restore its youthful state.

Biolife-Skinbooster also contains PDRN, a precious ingredient extracted from salmon roe, which has powerful repair and regenerative capabilities. It can stimulate the skin's own repair mechanism, promote cell regeneration, improve skin structure and enhance elasticity, thereby more effectively delaying skin ageing. The unique efficacy of PDRN can perfectly synergise with hyaluronic acid, bringing a comprehensive rejuvenating effect to the skin.




# Initial Treatment Stage



## Repair Damaged Dermal Layer

It is recommended to undergo 1 to 2 initial treatments, with a 2 to 4 week interval between treatments. The purpose of this stage is to repair the damaged dermal layer and improve the skin's firmness and elasticity.



## Immediate Tightening Effect

By injecting a high concentration of transparent hyaluronic acid gel, the skin can be effectively replenished with moisture, resulting in an immediate tightening effect.




## Continuous Maintenance and Improvement

Six months after the initial treatment, it is recommended to undergo a follow-up treatment. By injecting an appropriate amount of transparent hyaluronic acid gel, the skin can continue to be replenished with moisture, enhancing skin elasticity and reducing the formation of new wrinkles.

The Biolife-Skinbooster injection treatment can be customised based on individual needs. Our medical team will develop a targeted treatment plan, including the injection dosage and frequency, based on each person's skin condition, damage characteristics, and desired results, to ensure the optimal outcome of the treatment.

# The Unique Technology of Biolife-Skinbooster



## MDM Technology

Patented multi-crosslinking MDM technology using high-purity DVS.



## MBMT Technology

Monophasic microsphere technology, enhancing viscoelasticity.



## DMCT Technology

Divinyl sulfone multi-crosslinking technology, extending the lifetime of use.



## MDAP Technology

Multi-stage amphoteric purification, ensuring high purity and safety.

# **Biolife-Skinbooster - An Innovative Injection Treatment Certified with EU CE Mark**

Biolife-Skinbooster is an advanced injection treatment product that has obtained the EU CE certification. The CE mark is the strictest medical standard in Europe, indicating that Biolife-Skinbooster has reached the highest level of safety, efficacy and quality. This demonstrates its excellent medical standards and international status.

## **Widely Acclaimed**

Biolife-Skinbooster has been launched in multiple countries globally and has received widespread acclaim from professional doctors and patients. We will continue to rely on international standards to constantly improve product performance, providing better solutions for those seeking youthful and radiant skin.

# Product Specifications

Ingredients	Non-crosslinked high molecular weight hyaluronic acid 3,800m/kd + Low-crosslinked hyaluronic acid + PDRN
Volume	2 ml/syringe
Hyaluronic Acid Concentration	20mg/ml
Injection Strength	Viscosity 5-10N
Recommended Application Areas	Face, neck, neckline, hands, body
Needle	31g 13mm



# **Biolife-Skinbooster Restores Vitality to the Skin**

Dry, lacking elasticity skin regains its vitality under the care of Biolife-Skinbooster. Infusing nourishing nutrients, the moisturised, plump skin instantly radiates a dazzling glow. Fine lines and wrinkles gradually fade, replaced by a smooth, refined touch. Whether it's the face, neck or hands, Biolife-Skinbooster can infuse endless vitality into your skin.



# Why Choose Biolife-Skinbooster

## Highly Effective Hydration

Hyaluronic acid can absorb 1000 times its own weight in water, locking moisture in the skin to make it appear smoother and healthier.

## Safe and Reliable

Utilising patented technology and high-purity ingredients to ensure the safety and efficacy of the treatment.

## Long-Lasting Effects

Multi-crosslinking technology ensures the product's longevity in the skin, extending the treatment's benefits.

## Comprehensive Improvement

Not only enhances the skin's appearance, but also restores skin health from deep within, comprehensively improving skin quality.



# Production Process Comparison



## **BDDE Standard Production**

Each batch of 30,000 syringes requires 48 hours, including the preparation and dissolution of hyaluronic acid raw materials, crosslinking reaction using BDDE reagent, and strict temperature and pH control.



## **Washing and Filling**

BDDE production requires an additional 48 hours for thorough syringe washing and pre-filling and sealing of the product.



## **Biolife-Skinbooster Production**

Each batch of 2,000 syringes requires 168 hours, using a patented multi-crosslinking technology, with multiple purification and concentration processes to ensure product quality stability and safety.



## **Purification and Filling**

Biolife-Skinbooster production requires an additional 552 hours for multi-stage purification processes and aseptic pre-filling and sealing in a sterile environment, ensuring the safety and usability of each syringe.

# Biolife-Skinbooster

Biolife-Skinbooster provides a comprehensive and long-lasting solution, not only improving the appearance of the skin, but also restoring skin health from deep within, and comprehensively enhancing skin quality. Its unique manufacturing process and safe and reliable ingredients ensure the stable quality and safety of each syringe. Compared to regular hyaluronic acid injection products, Biolife-Skinbooster has invested more time and effort to ensure the exceptional quality of the final product.

## BIOLIFE SKIN BOOSTER





## **Exclusive Global Master Agent**

### **RICH COSMO (HK) LIMITED**

We are the exclusive global master agent for Biolife-Skinbooster, responsible for promoting and selling this innovative hyaluronic acid injection product. As a Hong Kong registered company, we have extensive market experience and a professional technical team, dedicated to providing high-quality services to customers worldwide.



## **Contact Us**

+852-39041158

[www.biolifeskinbooster.com](http://www.biolifeskinbooster.com)

[info@biolifeskinbooster.com](mailto:info@biolifeskinbooster.com)

Unit 802B, 8/F, Fook Heng Industrial Building, 1 Hok Yuen Street East, Hung Hom, Hong Kong